



BLOOD BANK REFRIGERATOR



HCBB-150L



Blood Bank Refrigerators are specialized medical devices that store blood and blood products at optimal temperatures to ensure their safety and quality for transfusions.

Glass Door with adhesive gasket

Protects samples and product temperature resistance.

High density polyurethane foam insulation

With an insulation thickness of $\pm 70\text{mm}$, it provides a good low (cold) value to maintain a cold temperature for a long time.

Automatic defrosting system with hot gas provides homogeneity and stability, accurate temperature precision.

Temperature Alarm System

High and low temperature alarm system for safe sample storage, door opening for temperature control in the room.

Non-CFC/HCFC cooling system. The compressor and refrigerant used are more environmentally friendly.

HCBB-1049L



Good Temperature



Good Insulation



Automatic defrosting system



Various alarm sistem



Non CFC

Alarm System

✓ Low Temperature

✓ High Temperature

✓ Door Open

MODEL		HCBB-150L	HCBB-1049L
(Approximately W x D x H mm)	Outer	$\pm 560 \times 720 \times 1415$	$\pm 1260 \times 890 \times 1975$
Weight (Approximately)		$\pm 105 \text{ kg}$	$\pm 250 \text{ kg}$
Capacity (Approximately)		150 L	1053 L
Shelf		3 EA	12 EA
Controller Display		7" full touchscreen	
Temperature	Range	2.0°C to 6.0°C	
Temperature Alarm		Low-High, Audible and Visual alarm	
Door		One door	Two door
Circulation		Fan	
Defrosting Technique		Automatic defrost with selenoid	
Safety		Over & Under Temp. Alarming, Thermal protector, Safety door lock	
Electric Supply		220 $\pm$ 10 VAC, 50/60 Hz	



PT. GOLDEN GLOBE MEDICA
providing a better and reliable solution



Head Office

Sentra Industri Terpadu Blok E2 No. 25, Pantai Indah Kapuk, Jakarta Utara 14470
Jakarta - Indonesia

Telp : +62 21 5698 2721 | +62 21 5698 3275
Fax : +62 21 5698 2720 | +62 21 5698 3276

Telp : +62 21 350 3011
Fax : +62 21 350 2667

Representative Office

• Yogyakarta
Sanggrahan No. 037 RT/RW.08/09 Banyuraden Gamping, Sleman
Yogyakarta 55293
Telp/ Fax : +62 274 530 4195